

Catalog	Description
HL - URCqPM2X	1.3ml

Store at **-20°C**  
 PI/HL-URCqPM2X - 02


## Product Description

**Rcore qPCR Mix (2X) with UDG** is a premix with concentrated reagents ready to use for Real-time Polymerase chain reactions (PCR). RNA template can be directly used for cDNA synthesis followed by amplification in a single tube which greatly reduces the total experimental time. The Rcore qPCR mix with UDG comes with all the necessary reagents that are required for cDNA synthesis and amplification in a single experimental setup. The 2X concentrated mix should be diluted to 1X with the template, primers, probes, and molecular grade water. In the regular PCR program, before the actual start of the amplification cycle, the inclusion of an additional step of 5 minutes at 53°C for cDNA synthesis facilitates the completion of the entire experiment in a single PCR Program.

The reagent along with Reverse transcriptase, HotStart Taq Polymerase enzyme, and dNTPs come with dUTPs and Uracil DNA Glycosylase (UDG) enzyme. The incorporation of dUTP in place of dTTP to the growing strand during PCR makes it an ideal degradation target for the UDG enzyme. The UDG enzyme gets inactivated at higher temperatures and therefore does not interfere with the amplification while clearing out any potential contamination that may be carried over from previous experiments. The mix is available with and without UDG.

## Features

- Contains Reverse transcriptase enzyme for cDNA synthesis, HotStart Taq Polymerase Enzyme, Uracil DNA Glycosylase (UDG) to provide a better amplification with 5'-3' Polymerase and Exonuclease Activity
- Includes buffer, and dNTPs for ready to use PCR applications
- High specificity and sensitivity
- Buffer enhancements guarantee performance and reliability



## Reagents supplied

Rcore qPCR Mix (2X) with UDG (1.3ml)

## Storage conditions

Rcore qPCR Mix (2X) with UDG (1.3ml) should be stored at **-20°C**

## Recommended reaction set-up for PCR

The Rcore PCR mix with UDG works the best for specific primers and is not recommended for random primers! Prepare PCRs using required volumes of freeze-thawed components in the PCR hood as recommended in the table below

### PCR Protocol

Component	Volume(μL)
Rcore q PCR Mix With UDG	13.0
Forward primer (10 pm/μL)	0.4-0.7
Reverse primer (10 pm/μL)	0.4-0.7
TaqMan Probe	Variable (user defined)
Template RNA or CDNA	Variable (user defined)
Nuclease Free Water	Variable (user defined)
<b>Total Volume</b>	<b>26</b>

The reaction setup is for guidance and it should be modified according to the user's needs. The current product is sufficient for 100 reactions if the above protocol is used

## Reference Guide

**Note:** Mix the vials with the reaction set up by gentle tapping. A short spin is recommended after gentle mixing to ensure that the reagent mix is not sticking to the walls of the PCR tube

### Thermocycling Protocol

Place the tubes in the Thermal Cycler and start the polymerase chain reaction protocol. Below is a general template for PCRs and should be optimized for good results

Recommended PCR Program			
Operation	Temp	Time	Cycles
For cDNA synthesis	53°C	5 min	1
Start of Amplification program			
Initial denaturation	95°C	15 min	1
	95°C	15 sec	35-40
Annealing & Extension	60°C	60 sec	cycles

### Applications

- End-point and quantitative RT-Polymerase chain reaction (PCR)
- Microarray analysis
- RACE PCR
- Amplification of RNA

#### Other PCR products from Huwel LifeSciences that you may be interested

S.No.	Product description	Catalogue No.
1.	Rcore qPCR Mix (2X)	HL - RCqPM2X
2.	Rcore qPCR Mix (4X)	HL - RCqPM4X
3.	Rcore qPCR Mix with UDG (4X) ***	HL-URCqPM4X
4.	Taq Polymerase Enzyme (5U/μl)	HL - Taq - 50 - 250Units HL - Taq - 100 - 500Units HL - Taq - 200 - 1000Units
5.	HotStart Taq Polymerase Enzyme (5U/μl)	HL - HSTaq - 50 - 250Units HL - HSTaq - 100 - 500Units HL - HSTaq - 200 - 1000Units
6.	QuickStart Taq Polymerase Enzyme (5U/μl)	HL - QSTaq - 50 - 250Units HL - QSTaq - 100 - 500Units HL - QSTaq - 200 - 1000Units

\*\*\* This reagent is available without Uracil DNA Glycosylase (UDG) also

For further information on protocols and details please contact our technical support:

[info@huwellifesciences.in](mailto:info@huwellifesciences.in)