

Quantiplus® Adenovirus Real-Time Quantitative PCR Kit



QT-ADV-25 : 25 rxns

QT-ADV-50 : 50 rxns QT-ADV-100 : 100 rxns

PI/QTADV-02

Intended Use

Quantiplus® Adenovirus Real-Time Quantitative PCR Kit is a Real-Time PCR based in vitro diagnostic assay for quantitation of Adenovirus in human plasma, CSF, and urine. The kit contains Amplification Mix with specific Primers and Probes, Standards (ADQS1-ADQS4) and Internal Control. The kit contains a second amplification system to identify possible PCR inhibition by using an exogenous internal control (Huwel IC-A Mix) without affecting the analytical sensitivity of the assay.

Background Information

Human Adenoviruses (HAdVs) are ubiquitous double-stranded DNA viruses that cause a wide array of diseases in humans including pharyngitis, pneumonia, gastroenteritis, hemorrhagic cystitis, and keratoconjunctivitis. They also cause life-threatening opportunistic infections in immunocompromised individuals and are responsible for outbreaks in certain populations. Diagnosis of HAdV is mostly done by cell culture or antigen detection methods.

Kit Components

Color Coding (Caps)	Contents	Description	25 rxns (QT-ADV-25)	50 rxns (QT-ADV-50)	100 rxns (QT-ADV-100)
Amber	Huwel Adeno Ready Mix	Probes and Primers for Adenovirus and Internal Control along with Amplification Mix	1 x 375μL	1 x 750μL	2 x 750μL
Natural	Huwel IC-A Mix	Internal Control	1 x 300 μL	1 x 600 μL	2 x 600 μL
Pink	Huwel ADQS1	2 X 10 ⁶ IU/μL	1 x 100 μL	1 x 100 μL	2 x 100 μL
Pink	Huwel ADQS2	2 X 10 ⁵ IU/μL	1 x 100 μL	1 x 100 μL	2 x 100 μL
Pink	Huwel ADQS3	2 X 10 ⁴ IU/μL	1 x 100 μL	1 x 100 μL	2 x 100 μL
Pink	Huwel ADQS4	2 X 10 ³ IU/μL	1 x 100 μL	1 x 100 μL	2 x 100 μL
White	Huwel PW	Purified Water	1 x 500 μL	1 x 500 μL	1 x 1 mL

Note: Please pay attention to the cap color coding and the tube contents.

Huwel PW: Molecular Biology Grade Purified Water.

Storage and Transportation Conditions

The kit should be transported at temperature below -20 °C. The kit is stable until the expiry date printed on the package, if the storage temperature is within -20 ± 5 °C. More than 4X freezing and thawing cycles reduces the assay sensitivity. For intermittent usage the reagents should be frozen in aliquots.

Technical specifications

Target Sequence	Conserved region in <i>Hexon</i> gene of Adenovirus
Specificity	100%
Sensitivity	0.5 IU/μL (250 IU/mL or 250 copies/mL)
Linear Range	$2X10^7 - 0.6 \text{ IU/}\mu\text{L} (1 \times 10^{10} - 3 \times 10^2 \text{IU/mL or } 1 \times 10^{10} - 3 \times 10^2 \text{ copies/mL})$
Reporting Units	IU/mL or Copies/ mL(1 IU = 1 copy)
Validated Specimen	Plasma
External Quality Assessment	QCMD EQA Panels



Assay Procedure

DNA Extraction

Quantiplus® Adenovirus Real-Time Quantitative PCR Kit has been validated using the following Viral DNA extraction kits: Recommended sample volume for extraction and elution are as follows:

S. No.	Name of the Extraction Kit	Recommended Sample volume for Extraction	Recommended Final Elution volume
1.	Huwel Nucleic Acid Extraction Kit - Version 2.0 (Cat. No. HL-NAX-100)	200 μL	100 μL
2.	QIAamp [®] DNA Blood Mini Kit (Cat. No. 51104)	200 μL	100 μL

Note: Customer can also validate their own extraction process using other Viral DNA extraction Kits. IC-A mix can be added at the extraction step or while setting up the PCR

qPCR Protocol Flow 1.5 mL Centrifuge tube PCR tube PCR tube PCR tube Real-Time PCR Result analysis # If IC not added during extraction step

Preparation of Reaction Master mix

Components	Volume per reaction (for 26μL)
Huwel Adeno Ready Mix	15.0
Huwel IC-A Mix (if not added at extraction step)	1.0
Extracted DNA/Huwel ADQS1- Huwel ADQS4/ Huwel PW	10.0

It is necessary to keep all components at +2 °C to +8 °C during the PCR preparation. Close the tubes and centrifuge briefly before proceeding to the thermal cycler.

Cycling Conditions

Steps	No. of cycles	Temperature (°C)	Time
1 (Initial denaturation)	1	95	15 min.
2 (PCR cycling)	45	95	15 sec.
2 (FCK Cycling)	43	60*	1 min

*Plate read/Data acquisition in FAM and YAKIMA YELLOW/ VIC/ HEX channels in Bio-Rad™, CFX 96 and Thermo QS5 Real-Time PCR System.



Sample analysis and Interpretation

The criteria for the acceptance of the assay should be met before the interpretation of the unknown sample results as described in Table 1 below and also the slope of the standard curve (standards in FAM channel) is between -3.1 to -3.6, (at least three standards should be included,) and PCR efficiency is between 90% to 110% (0.9 to 1.1). Interpret the results of unknown samples as mentioned in Table 2

Table 1:

Control	FAM (ADV)	YAKIMA YELLOW/ VIC/ HEX (IC)	
If Internal Control (IC-A Mix) is added during extraction			
Standards (ADQS1 to ADQS4)	٧	-	
Negative Control (NC)	-	-	
If Internal Control (IC-A Mix) is added during preparation of reaction master mix			
Standards (ADQS1 to ADQS4)	٧	V	
Negative Control (NC)	-	V	

Table 2:

S.No	FAM	YAKIMA YELLOW/	← Fluorophore Interpretation	Conclusion
	(ADV)	VIC/HEX (IC)	interpretation	
1	٧	٧	Adeno DNA detected within quantitation	
2	٧	-	range	Proceed for further Analysis
3	-	٧	Adeno DNA below quantitation limit	
4	_	_	Possible inhibition of PCR	Dilute the DNA sample (1:10)
4	- Fossible illimbition of FCK	and repeat the Assay		

Viral load calculation (Conversion of IU/μL to IU/mL)

IU/mL = <u>Obtained IU/μL X Elution Volume (μL)</u>
Sample volume in mL

For calculating the result of diluted sample (1:10); multiply the observed IU/mL value by dilution factor, 10

Result (IU/μL) x Elution Volume (μL)
Result of 1:10 diluted sample (IU/mL) = Dilution Factor (10) x ------Sample Volume (mL)



Reporting Comments

Results in IU/mL	Comments	
Target not detected	Adeno DNA not detected in the given sample	
	Adeno DNA detected but below the lower limit of the	
	linear range of the assay. The reproducibility of the positive	
<300	result is not assured	
300 to 1X 10 ¹⁰	Adeno DNA detected within the linear range of the assay	
>1 X 10 ¹⁰	Adeno DNA detected but above linear range of the assay,	
	dilute the sample and repeat the assay for accurate result	

Validated Instruments

- Thermo QS5 Real-Time PCR System
- Bio-Rad[™] CFX 96



HLSS Manufacturing Pvt Ltd Plot No's M14, M15, M16, TSIIC Medical device park Sultanpur village, Ameenpur Mandal, SangareddyDist, TS-502319