

## Quantiplus® HCV FAST RT PCR Kit (Real-Time Quantitative PCR Kit)



QTF-HCV-25 : 25 rxns  
 QTF-HCV-50 : 50 rxns  
 QTF-HCV-100 : 100 rxns



PI/QTFHCV-07

### Intended Use

Quantiplus® HCV FAST RT-PCR Kit is a reverse transcription Real-Time PCR based in vitro diagnostic assay for quantitation of Hepatitis C Virus (1-6 genotypes) in human plasma. The kit contains single tube qPCR mix with Reverse Transcriptase and UDG/UNG, Primer Probe Mix (PPM), Standards (HCVFQS1-HCVFQS4), and Internal Control (IC-B mix). This advanced formulation enables performance of fast PCR in shorter run time ( $\leq 60$  min), and UDG/UNG helps in controlling PCR carryover contamination.

**This kit is not to be used for screening of blood/blood products from blood donors.**

### Background Information

Hepatitis C virus causes liver disease which can lead to both acute and chronic hepatitis. It can cause mild illness lasting a few weeks to a serious, lifelong illness. A significant number of chronically infected individuals develop cirrhosis or liver cancer. The test is used for HCV Viral load quantitation as it is important to start/ for follow up therapy.

### Kit Components

Color Coding (Caps)	Contents	Description	25 rxns (QTF-HCV-25)	50 rxns (QTF-HCV-50)	100 rxns (QTF-HCV-100)
Green	RNA Fast qPCR Mix with UDG/UNG (4X)	PCR Amplification mix	1 x 165 $\mu$ L	1 x 330 $\mu$ L	2 x 330 $\mu$ L
Amber	HCV Fast PPM	Target specific Primer Probe Mix	1 x 50 $\mu$ L	1 x 100 $\mu$ L	2 x 100 $\mu$ L
Natural	IC-B Mix	Exogenous Internal Control-B mix	1 x 300 $\mu$ L	1 x 600 $\mu$ L	2 x 600 $\mu$ L
Pink	HCVFQS1	1 X 10 <sup>4</sup> IU/ $\mu$ L	1 x 100 $\mu$ L	1 x 100 $\mu$ L	2 x 100 $\mu$ L
Pink	HCVFQS2	1 X 10 <sup>3</sup> IU/ $\mu$ L	1 x 100 $\mu$ L	1 x 100 $\mu$ L	2 x 100 $\mu$ L
Pink	HCVFQS3	1 X 10 <sup>2</sup> IU/ $\mu$ L	1 x 100 $\mu$ L	1 x 100 $\mu$ L	2 x 100 $\mu$ L
Pink	HCVFQS4	1 X 10 <sup>1</sup> IU/ $\mu$ L	1 x 100 $\mu$ L	1 x 100 $\mu$ L	2 x 100 $\mu$ L
White	MBGPW	Purified water	2 x 500 $\mu$ L	2 x 1 mL	2 x 1 mL

**Note: Please pay attention to the cap color coding and the tube contents.**

MBGPW: Molecular Biology Grade Purified Water

### Storage and Transportation Conditions

The kits should be transported at temperature below  $-20^{\circ}\text{C}$ . The kit is stable until the expiry date printed on the package, if the storage temperature is within  $-20 \pm 5^{\circ}\text{C}$ . More than 4X freezing and thawing cycles reduces the assay sensitivity. For intermittent usage the reagents should be frozen in aliquots.

### Technical specifications

Target Sequence	Conserved region of 5' UTR
Specificity	1-6 genotypes with 100% specificity
Sensitivity	$6.0 \times 10^{-2}$ IU/ $\mu$ L ( $3.0 \times 10^1$ IU/mL or $6.5 \times 10^1$ copies/mL)
Linear Range	7.3 log IU/mL ( $2.25 \times 10^7$ IU/ $\mu$ L) to 1.68 log/mL ( $9.6 \times 10^{-2}$ IU/ $\mu$ L)

Reporting Units	IU/μL (1 IU = 2.18 copies)
Validated Specimen	Plasma
External Quality Assessment	QCMD EQA Panels
Standards	Calibrated against to HCV NIBSC code: 06/102


**Assay Procedure**

**RNA Extraction**

Quantiplus® HCV FAST RT PCR Kit (Real-Time Quantitative PCR Kit) has been validated using the Viral RNA extraction kits mentioned below. Recommended sample volume for extraction and elution are as follows:

S. No.	Name of the Extraction Kit	Recommended Sample volume for Extraction	Recommended Final Elution volume
1.	Huwel Nucleic Acid Extraction Kit - Version 2.0 (Cat. No. HL-NAX-100)	200 μL	100 μL
2.	QIAamp® Viral RNA Mini Kit (Cat. No. 52904)	140 μL	60 μL

*Note: Customer can also validate their own extraction process using other Viral RNA extraction Kits. IC-B mix can be added at the extraction step or while setting up the PCR*



**Preparation of Reaction Master mix**

Components	Volume per reaction (for 25μL)
RNA Fast qPCR Mix with UDG/UNG (4X)	6.5
HCV Fast PPM	2.0
IC-B Mix (if not added at extraction step)	1.0
MBGPW	5.5
Extracted RNA/ HCVFQS1- HCVFQS4 / MBGPW	10.0

It is necessary to keep all components at +2 °C to +8 °C during the PCR preparation. Close the tubes and centrifuge briefly before proceeding to the thermal cycler.

**Cycling Conditions**

Steps	No. of cycles	Temperature (°C)	Time
1(Reverse Transcription)	1	53	5 min.
2 (Initial denaturation)	1	95	1 min.
3 (PCR cycling)	45	95	10 sec.
		60*	10 sec.

\*Plate read/Data acquisition in **FAM** and **TEXAS RED** channels in Bio-Rad™ CFX 96. For Thermo Q55 Real-Time PCR System, use **FAM** and **ROX** channels. For Rotor-Gene Q 5 plex, use **Green** and **Orange** channels.

**Sample analysis and Interpretation**

The criteria for the acceptance of the assay should be met before the interpretation of the unknown sample results as described in Table 1 below and also ensure that the slope of the standard curve (standards in FAM channel) is between -3.1 to -3.6 (at least three standards should be included) and PCR efficiency is between 90% to 110% (0.9 to 1.1). Interpret the results of unknown samples as mentioned in Table 2

**Table 1:**

Control	FAM (HCV)	TEXAS RED (IC)
If Internal Control (IC-B Mix) is added during extraction		
Standards (HCVFQS1 to HCVFQS4)	√	-
Negative Control (NC)	-	-
If Internal Control (IC-B Mix) is added during preparation of reaction master mix		
Standards (HCVFQS1 to HCVFQS4)	√	√
Negative Control (NC)	-	√

**Table 2:**

S.No	FAM (HCV)	TEXAS RED (IC)	Interpretation	Conclusion
1	√	√	HCV RNA detected within quantitation range	Proceed for further Analysis
2	√	-		
3	-	√	HCV RNA below quantitation limit	Dilute the RNA sample (1:10) and repeat the Assay
4	-	-	Possible inhibition of PCR	

**Viral load calculation (Conversion of IU/μL to IU/mL)**

$$\text{IU/mL} = \frac{\text{Obtained IU/}\mu\text{L} \times \text{Elution Volume } (\mu\text{L})}{\text{Sample volume in mL}}$$

For calculating the result of diluted sample (1:10); multiply the observed IU/mL value by dilution factor, 10

$$\text{Result of 1:10 diluted sample (IU/mL)} = \text{Dilution Factor (10)} \times \frac{\text{Result (IU/}\mu\text{L)} \times \text{Elution Volume } (\mu\text{L})}{\text{Sample Volume (mL)}}$$

**Reporting comments**

Results in IU/mL	Comments
Target not detected	HCV RNA not detected in the given sample
<48	HCV RNA detected but below the lower limit of the quantitation range of the assay. The reproducibility of the positive result is not assured
48 to $1.1 \times 10^{10}$	HCV RNA detected within the linear range of the assay
$\geq 1.1 \times 10^{10}$	HCV RNA detected but above linear range of the assay, dilute the sample for accurate result.

**Validated Instruments**

- Thermo QS5 Real-Time PCR System
- Bio-Rad™ CFX 96
- Rotor-Gene Q 5 plex



HLSS Manufacturing Pvt Ltd  
 Plot No's M14, M15, M16, TSIIIC Medical device park  
 Sultanpur village, Ameenpur Mandal,  
 SangareddyDist, TS-502319