

Huwel Genomic DNA Extraction Kit



HL-GDX-100: 100 Exts

PI/HLGDX-00

Intended use

Huwel Genomic DNA Extraction Kits are intended for molecular biology applications such as PCR, RT-PCR, and sequencing. These kits are not intended for the diagnosis, prevention, or treatment of a disease.

Background Information

The Huwel Genomic DNA Extraction Kit is a fast, economical and easy method of genomic DNA isolation from blood or tissue. The buffer system provided in the kit allows efficient removal of RBC and lysis of WBC followed by selective binding of nucleic acids to the spin column. These kits produce high quality and yield of genomic DNA.

Kit components

Component	Quantity	Storage
10X RBC Lysis Buffer	50 mL	RT
Spin Column	100	RT
Collection tube	200	RT
Huwel Lysis Buffer	40 mL	RT
Huwel Wash Buffer 1	50 mL	RT
Huwel Wash Buffer 2	2 x 50 mL	RT
Huwel Elution Buffer	20 mL	RT
Lyophilized Poly A	500 µg	-15°C to -25°C
Proteinase K	100 mg	-15°C to -25°C
Kit Protocol	1	NA

Note: Huwel Wash Buffer 1 and Huwel Wash Buffer 2 are provided with Ethanol. No need to add additional Ethanol.

Add 10 µl β-ME or 20 µl of 2 M dithiothreitol (DTT) per 1 ml of Huwel Lysis Buffer before use.

Storage

- Huwel Genomic DNA Extraction Kits are stable at room temperature (15–25°C) until the expiration date mentioned on the kit
- **Store the Lyophilized Poly A at 4°C and Lyophilized Proteinase K at -15°C to -25°C upon receipt of the kit**
- Reconstitute the Lyophilized Poly A in 500 µL of Huwel Elution Buffer. Prepare 50 µL aliquots in sterile vials and store at -15°C to -25°C; stable for 12 months
- Reconstitute the Lyophilized Proteinase K in 5 mL of Huwel Elution Buffer. Prepare small aliquots according to the user requirement and store at -15°C to -25°C; stable for 12 months

Huwel Genomic DNA Extraction Kit

Procedure

RBC lysis from whole blood:

- Mix 1 volume of whole blood with 5 volumes of 1x RBC Lysis Buffer.
- Incubate for 10–15 min on ice. Mix by vortexing briefly 2 times during incubation.
- Centrifuge at 400 x g for 10 min at 4°C, and completely remove and discard supernatant.
- Add 2 volumes of 1x RBC Lysis Buffer to the cell pellet. Resuspend cells by vortexing briefly.
- Centrifuge at 400 x g for 10 min at 4°C, and completely remove and discard supernatant.
- Add 400 µl of Huwel Lysis buffer to the cell pellet. Mix thoroughly by vortex or pipetting to remove any clumps and proceed to Nucleic Acid extraction.

RBC lysis from tissue or solid tumor samples:

- Dissociate tissue or solid tumor samples into single cells.
- Centrifuge cells at 400 x g for 5 minutes at room temperature and remove the supernatant carefully.
- Resuspend cell pellet in 5-10 mL of 1X RBC Lysis Buffer.
- Incubate for 10 minutes at room temperature.
- Centrifuge at 400 x g for 5 minutes at room temperature and remove the supernatant
- carefully.
- Re-suspend the cell pellet in appropriate volume of 1x PBS. Proceed to Nucleic Acid extraction.

Huwel Genomic DNA Extraction Kit

Nucleic Acid extraction:

1. To 400 μ L of cell lysate in Huwel Lysis Buffer, add 5 μ L of **Poly A** and 50 μ L of **Proteinase K** in a microfuge tube.
2. Mix well and incubate for 10 min at 72°C
3. Add 100 μ L of **Huwel Lysis Buffer**, vortex the tube for 5 sec and briefly spin to collect the sample droplets inside the cap
4. Load the entire sample lysate (555 μ L) into Huwel spin column-collection tube assembly
5. Centrifuge at 8,000 x g for 1 min. Discard the flowthrough.
6. Change the collection tube, add 500 μ L of **Huwel Wash Buffer 1**
7. Centrifuge 8,000 x g for 1 min. Discard the flowthrough
8. Add 450 μ L of **Huwel Wash Buffer 2**
9. Centrifuge 8,000 x g for 1 min and discard the flowthrough
10. Add 450 μ L of **Huwel Wash Buffer 2**
11. Centrifuge 8,000 x g for 1 min. Discard the flowthrough
12. Change the collection tube (microfuge tube not provided), give empty spin by centrifuging at 13,000 x g for 1 min
13. Place the spin column into a fresh microfuge tube
14. Add 50/100 μ L of **Huwel Elution Buffer** into the spin column
15. Incubate the spin column at RT for 5 mins
16. Centrifuge 8,000 x g for 1 min
17. Store the eluent at -20°C

Huwel Genomic DNA Extraction Kit

Troubleshooting guide

Problem	Possible cause	Recommendation
Low nucleic acid quality or yield	Buffers stored at sub-optimal conditions	Store the kit at room temperature (15–25°C). Close reagent bottles tightly to avoid contamination and maintain stability. Store the reconstituted Poly A and Proteinase K at -15°C to -25°C.
	Incomplete mixing of sample and reagents	Mix the sample tube well after adding each reagent.
Poor elution of Nucleic acids	Drift in pH of water	If an alternate source of water or buffer is used for elution, make sure that the pH matches with the Elution buffer provided in the kit.
Low RNA yield	RNase contamination	Work in RNase-free environment. Use fresh samples or store the samples at -80°C until they can be processed. Use eluted RNA immediately for downstream applications or store at -80°C.



HLSS Manufacturing Pvt Ltd
 Plot No's M14, M15, M16, TSIC Medical device park
 Sultanpur villiage, Ameenpur Mandal,
 Sangareddy Dist, TS-502319